

Anti-Fndc5 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Nau1

Item Number: CM0023092



Introduction

Anti-FNDc5 rabbit polyclonal antibody validated for WB, ELISA. Reacts with human, rat. Targets fibronectin type III domain-containing protein 5 (irisin), a key exercise-induced hormone.

[Learn More](#)

Parameter	Value
Product Name	FNDc5 Rabbit pAb
Remarks/Alias	FRCP2; irisin; FNDc5
Species	Human
GeneID (Human)	252995
GeneID	252995
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 113-212 of human FNDc5 (NP_715637.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q8NAU1
Protein Weight	24kDa
Shipping	Ice bag