

Anti-Fnta Polyclonal Antibody For Wb, Elisa - P49354

Item Number: CM0024724



Introduction

Rabbit polyclonal antibody against human FNTA, validated for WB and ELISA. Cross-reacts with mouse and rat. Target involved in protein prenylation, relevant for cancer and neuromuscular research.

[Learn More](#)

Parameter	Value
Product Name	FNTA Rabbit pAb
Remarks/Alias	FPTA, PTAR2, PGGT1A, FNTA
Species	Human
GeneID (Human)	2339
GeneID	2339
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 280-379 of human FNTA (NP_002018.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49354
Protein Weight	44kDa
Shipping	Ice bag