

## Anti-Folate Binding Protein (Folr1) Rabbit Polyclonal Antibody For Wb - P15328

Item Number: CM0023521



### Introduction

Rabbit polyclonal antibody against human FOLR1 (folate receptor alpha) for WB, IHC-P, IF/ICC, ELISA. Validated in human, mouse, rat. Ideal for cancer and folate metabolism studies.

[Learn More](#)

| Parameter        | Value                                                          |
|------------------|----------------------------------------------------------------|
| Product Name     | Folate Binding Protein(FBP) / FOLR1 Rabbit pAb                 |
| Remarks/Alias    | FBP; FOLR; NCFTD; FRalpha; Folate Binding Protein(FBP) / FOLR1 |
| Species          | Human                                                          |
| Gene ID (Human)  | 2348                                                           |
| Gene ID          | 2348                                                           |
| Immunogen        | Synthetic peptide                                              |
| Source           | Rabbit                                                         |
| Category         | Polyclonal Antibodies                                          |
| Application      | WB, IHC-P, IF/ICC, ELISA                                       |
| Cross Reactivity | Human, Mouse, Rat                                              |
| SWISS            | P15328                                                         |
| Protein Weight   | 30kDa                                                          |
| Shipping         | Ice bag                                                        |