

Anti-Fopnl Polyclonal Antibody For Wb, Elisa - Q96Nb1

Item Number: CM0024205

Introduction

High-quality rabbit polyclonal antibody targeting FOPNL (Centrosomal protein 20) for Western blot and ELISA applications. Validated for human and mouse samples, ideal for cilia and centrosome research.

[Learn More](#)

Parameter	Value
Product Name	FOPNL Rabbit pAb
Remarks/Alias	FOPNL; FOR20; C16orf63; PHSECRG2
Species	Human
GeneID (Human)	123811
GeneID	123811
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-100 of human FOPNL (NP_653201.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96NB1
Protein Weight	20kDa
Shipping	Ice bag