


Anti-Fos B Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P53539

Item Number: CM0007703



Introduction

High-quality anti-Fos B rabbit monoclonal antibody validated for Western blot, IHC-P, and ELISA. Recognizes human Fos B (AP-1 transcription factor) with a molecular weight of 36kDa. Ideal for studies of transcriptional regulation, T cell biology, and neurobiology.

[Learn More](#)

Parameter	Value
Product Name	Anti-Fos B Rabbit Monoclonal Antibody
Remarks/Alias	AP-1; G0S3; G0S3; G0SB; Fos B
Species	Human
GeneID (Human)	2354
GeneID	2354
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Fos B (NP_006723.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P53539
Protein Weight	36kDa
Shipping	Ice bag