


Anti-Fos Rabbit Polyclonal Antibody For Wb, Elisa - P01100

Item Number: CM0024104



Introduction

Rabbit polyclonal antibody targeting c-Fos (AP-1 transcription factor subunit). Reacts with human, mouse, rat. Suitable for WB and ELISA; ideal for signal transduction, proliferation, and differentiation research.

[Learn More](#)

Parameter	Value
Product Name	FOS Rabbit pAb
Remarks/Alias	p55; AP-1; C-FOS; FOS
Species	Human, Mouse, Rat
GeneID (Human)	2353
GeneID	2353
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FOS (NP_005243.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01100
Protein Weight	41kDa
Shipping	Ice bag