


Anti-Fosb Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P53539

Item Number: CM0016835



Introduction

High-quality rabbit polyclonal antibody against human FOSB, validated for WB, IHC-P, ELISA. Ideal for studying AP-1 transcription factor functions in nucleosomal enhancer regulation and neurogenesis.

[Learn More](#)

Parameter	Value
Product Name	FOSB Rabbit pAb
Remarks/Alias	AP-1; G0S3; GOS3; GOSB; FOSB
Species	Human
GeneID (Human)	2354
GeneID	2354
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FOSB (NP_006723.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53539
Protein Weight	36kDa
Shipping	Ice bag