

Anti-Fosl1 (Kd Validated) Rabbit Polyclonal Antibody For Wb, Elisa - P15407

Item Number: CM0024682

Introduction

KD validated rabbit polyclonal antibody against human FOSL1 (FRA1) for Western blot and ELISA applications. Sourced from CamelBio's one-stop IVD raw material platform, this antibody detects a ~29 kDa band. Swiss-Prot: P15407.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] FOSL1 Rabbit pAb
Remarks/Alias	FRA; FRA1; fra-1; [KD Validated] FOSL1
Species	Human
GeneID (Human)	8061
GeneID	8061
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human FOSL1 (NP_005429.1?).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P15407
Protein Weight	29kDa
Shipping	Ice bag