

Anti-Fosl2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P15408

Item Number: CM0030213

Introduction

Rabbit polyclonal antibody against human FOSL2 (FRA2), validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human and rat samples. Suitable for studying transcription regulation and osteoclast biology. - P15408

[Learn More](#)

Parameter	Value
Product Name	FOSL2 Rabbit pAb
Remarks/Alias	FRA2; FOSL2
Species	Human
GeneID (Human)	2355
GeneID	2355
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 187-326 of human FOSL2 (NP_005244.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	P15408
Protein Weight	35kDa
Shipping	Ice bag