

Anti-Foxa1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P55317

Item Number: CM0024426



Introduction

Rabbit polyclonal antibody targeting human FOXA1 (HNF3A/TCF3A), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for studying pioneer factor function and hormone receptor signaling in development and cancer research.

[Learn More](#)

Parameter	Value
Product Name	FOXA1 Rabbit pAb
Remarks/Alias	HNF3A, TCF3A, FOXA1
Species	Human
GeneID (human)	3169
GeneID	3169
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 263-472 of human FOXA1 (NP_004487.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55317
Protein Weight	49kDa
Shipping	Ice bag