

Anti-Foxc1 Polyclonal Antibody For Wb, Elisa - Q12948

Item Number: CM0025947

Introduction

FOXC1 Rabbit Polyclonal Antibody validated for Western blot and ELISA, reactive with human, mouse, and rat. Suitable for studying transcription factor function in development, cancer, and ocular research.

[Learn More](#)

Parameter	Value
Product Name	FOXC1 Rabbit pAb
Remarks/Alias	ARA; IGDA; IHG1; ASGD3; FKHL7; IRID1; RIEG3; FREAC3; FREAC-3; FOXC1
Species	Human
GeneID (human)	2296
GeneID	2296
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 404-553 of human FOXC1 (NP_001444.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12948
Protein Weight	57kDa
Shipping	Ice bag