

Anti-Foxd1 Rabbit Polyclonal Antibody For Wb, Elisa - Q16676

Item Number: CM0024554

Introduction

Rabbit polyclonal antibody targeting human FOXD1, validated for Western blot and ELISA. Cross-reacts with mouse and rat. Suitable for transcription factor studies.

[Learn More](#)

Parameter	Value
Product Name	FOXD1 Rabbit pAb
Remarks/Alias	FKHL8; FREAC4; FREAC-4; FOXD1
Species	Human
Gene ID (Human)	2297
Gene ID	2297
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FOXD1 (NP_004463.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16676
Protein Weight	46kDa
Shipping	Ice bag