

Anti-Foxi1 Polyclonal Antibody For Wb, Elisa - Q12951

Item Number: CM0020245

Introduction

Polyclonal antibody targeting FOXI1, a forkhead transcription factor essential for inner ear and kidney development. Validated for WB and ELISA in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	FOXI1 Rabbit pAb
Remarks/Alias	HFH3, FKHL10, HFH-3, FKHL10, FREAC6, FREAC-6, FOXI1
Species	Human
GeneID (Human)	2299
GeneID	2299
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human FOXI1 (NP_036320.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q12951
Protein Weight	41kDa
Shipping	Ice bag