

Anti-Foxj1 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q92949

Item Number: CM0020071

Introduction

Rabbit polyclonal antibody targeting human FOXJ1, a transcription factor essential for motile cilia formation. Validated for WB, IP, ELISA applications, with cross-reactivity to mouse and rat. Ideal for ciliopathy and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	FOXJ1 Rabbit pAb
Remarks/Alias	HFH4; HFH-4; CILD43; FKHL13; FOXJ1
Species	Human
GeneID (human)	2302
GeneID	2302
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-421 of human FOXJ1. (NP_001445.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q92949
Protein Weight	45kDa
Shipping	Ice bag