

Anti-Foxk1 Rabbit Polyclonal Antibody For Wb, Elisa - P85037

Item Number: CM0014534



Introduction

Rabbit polyclonal antibody against human FOXK1, validated for WB and ELISA. FOXK1 is a forkhead transcription factor regulating glucose metabolism, muscle differentiation, and autophagy. Suitable for research applications.

[Learn More](#)

Parameter	Value
Product Name	FOXK1 Rabbit pAb
Remarks/Alias	FOXK1L; FOXK1
Species	Human
Gene ID (Human)	221937
Gene ID	221937
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 600-733 of human FOXK1 (NP_001032242.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P85037
Protein Weight	75kDa
Shipping	Ice bag