

Anti-Foxk2 Rabbit Polyclonal Antibody For Wb And Elisa - Q01167

Item Number: CM0016682



Introduction

FOXK2 Rabbit Polyclonal Antibody (pAb) for WB and ELISA. Detects human, mouse, and rat FOXK2. Ideal for transcription regulation, metabolism, and autophagy research. SWISS: Q01167.

[Learn More](#)

Parameter	Value
Product Name	FOXK2 Rabbit pAb
Remarks/Alias	ILF; ILF1; ILF-1; nGTBP; FOXK2
Species	Human
GeneID (Human)	3607
GeneID	3607
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 360-430 of human FOXK2 (NP_004505.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01167
Protein Weight	69kDa
Shipping	Ice bag