

Anti-Foxl2 Monoclonal Antibody For Wb, Elisa - P58012

Item Number: CM0008289



Introduction

Rabbit monoclonal FOXL2 antibody for Western blot and ELISA. Cross-reacts with human, mouse, rat FOXL2. Suitable for research in ovarian differentiation, transcription regulation, and disease mechanisms.

[Learn More](#)

Parameter	Value
Product Name	FOXL2 Rabbit mAb
Remarks/Alias	BPES; PFRK; POF3; BPES1; PINTO; FOXL2
Species	Human
GeneID (Human)	668
GeneID	668
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FOXL2 (P58012)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P58012
Protein Weight	39kDa
Shipping	Ice bag