

Anti-Foxm1 Polyclonal Antibody For Wb, If-P, Elisa - Q08050

Item Number: CM0019755



Introduction

High-quality rabbit polyclonal antibody against FOXM1, validated for WB, IF-P, and ELISA. Detects human, mouse, and rat FOXM1. Recombinant immunogen corresponding to aa 29-128 of human FOXM1. Ideal for cell cycle, DNA repair, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	FOXM1 Rabbit pAb
Remarks/Alias	MPP2; HFH11; HNF-3; INS-1; MPP-2; PIG29; FKHL16; FOXM1A; FOXM1B; FOXM1C; HFH-11; TRIDENT; MPHOSPH2; FOXM1
Species	Human
GeneID (Human)	2305
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-128 of human FOXM1 (NP_068772.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q08050
Protein Weight	84kDa
Shipping	Ice bag