

Anti-Foxn2 Rabbit Monoclonal Antibody For Wb, Elisa - P32314

Item Number: CM0006404

Introduction

Rabbit monoclonal antibody against human FOXN2 (Forkhead box protein N2) for WB and ELISA applications. Validated by knockout, recognizes human FOXN2. Suitable for HTLV-I LTR binding studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] FOXN2 Rabbit mAb
Remarks/Alias	HTLF; N2
Species	Human
GeneID (human)	3344
GeneID	3344
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 204-431 of human FOXN2. (NP_002149.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P32314
Protein Weight	47kDa
Shipping	Ice bag