

Anti-Foxn2 Rabbit Polyclonal Antibody For If/Icc, Elisa - P32314

Item Number: CM0029819

Introduction

FOXN2 Rabbit Polyclonal Antibody validated for IF/ICC and ELISA. Reacts with human FOXN2, a transcription factor binding HTLV-I LTR. Ideal for studies of viral gene regulation, transcription control, and FOX family biology.

[Learn More](#)

Parameter	Value
Product Name	FOXN2 Rabbit pAb
Remarks/Alias	HTLF; FOXN2
Species	Human
GeneID (Human)	3344
GeneID	3344
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 172-431 of human FOXN2 (NP_002149.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P32314
Protein Weight	47kDa
Shipping	Ice bag