

Anti-Foxn2 Rabbit Polyclonal Antibody For Wb And Elisa - P32314

Item Number: CM0028312

Introduction

FOXN2 rabbit polyclonal antibody validated for Western blotting and ELISA. Reacts with human and mouse FOXN2 (HTLF), a transcription factor that binds HTLV-I LTR. Suitable for chromatin and viral transcription research.

[Learn More](#)

Parameter	Value
Product Name	FOXN2 Rabbit pAb
Remarks/Alias	HTLF; FOXN2
Species	Human
Gene ID	3344
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 172-431 of human FOXN2 (NP_002149.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P32314
Protein Weight	47kDa
Shipping	Ice bag