

## Anti-Foxn3 Polyclonal Antibody For Wb And Elisa - 000409

Item Number: CM0015177



### Introduction

High-quality rabbit polyclonal antibody against FOXN3, a transcriptional repressor involved in DNA damage checkpoints. Validated for WB and ELISA, with cross-reactivity to mouse. Ideal for cell cycle and DNA damage research.

[Learn More](#)

| Parameter        | Value                                                                                                                                   |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | FOXN3 Rabbit pAb                                                                                                                        |
| Remarks/Alias    | CHE31; PRO1635; C14orf116; FOXN3                                                                                                        |
| Species          | Human                                                                                                                                   |
| GeneID (Human)   | 1112                                                                                                                                    |
| GeneID           | 1112                                                                                                                                    |
| Immunogen        | Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 269-468 of human FOXN3 (NP_005188.2). |
| Source           | Rabbit                                                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                                                   |
| Application      | WB, ELISA                                                                                                                               |
| Cross Reactivity | Mouse                                                                                                                                   |
| SWISS            | O00409                                                                                                                                  |
| Protein Weight   | 54kDa                                                                                                                                   |
| Shipping         | Ice bag                                                                                                                                 |