

## Anti-Foxo1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q12778

Item Number: CM0030287



### Introduction

Rabbit polyclonal antibody against FOXO1 (FKHR), validated for Western blot, immunofluorescence/immunocytochemistry, and ELISA in human, mouse, rat, and monkey samples. Ideal for diabetes, metabolism, and apoptosis research.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | FOXO1 Rabbit pAb                                                                                                    |
| Remarks/Alias    | FKH1, FKHR, FOXO1A, FOXO1                                                                                           |
| Species          | Human                                                                                                               |
| GeneID (Human)   | 2308                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 340-580 of human FOXO1 (NP_002006.2). |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, IF/ICC, ELISA                                                                                                   |
| Cross Reactivity | Human, Mouse, Rat, Monkey                                                                                           |
| SWISS            | Q12778                                                                                                              |
| Protein Weight   | 70kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |