

Anti-Foxo1/Foxo3A/Foxo4 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q12778 / O43524 / P98177

Item Number: CM0020869



Introduction

Rabbit polyclonal antibody recognizing FOXO1, FOXO3A, and FOXO4; suitable for WB, IHC-P, and ELISA with human, mouse, and rat samples. Ideal for metabolism, apoptosis, and transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	FOXO1/FOXO3A/FOXO4 Rabbit pAb
Species	Human
GeneID (Human)	2308/2309/4303
GeneID	2308/2309/4303
Immunogen	Synthetic peptide; A synthetic peptide of human FOXO1/FOXO3A/FOXO4 (NP_002006.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12778/O43524/P98177
Shipping	Ice bag