

Anti-Foxo3A Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - Q9Wvh4

Item Number: CM0025304

Introduction

FOXO3A Rabbit Polyclonal Antibody detects mouse and rat FOXO3A. Suitable for WB, IF/ICC, IF-P, and ELISA applications in apoptosis, autophagy, and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	FOXO3A Rabbit pAb
Remarks/Alias	Fkhr2; FKHL1; Foxo3a; 1110048B16Rik; 2010203A17Rik
Species	Mouse
GeneID (Human)	2309
GeneID	56484
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 371-470 of mouse FOXO3A (NP_062714.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9WVH4
Protein Weight	71kDa
Shipping	Ice bag