

Anti-Foxo3A Rabbit Monoclonal Antibody For Wb, Elisa - O43524

Item Number: CM0012000



Introduction

Rabbit monoclonal antibody targeting FOXO3A, validated for WB and ELISA in human and mouse samples. Recognizes a 71 kDa protein; detects endogenous levels. Ideal for apoptosis, autophagy, and oxidative stress research.

[Learn More](#)

Parameter	Value
Product Name	FOXO3A Rabbit mAb
Remarks/Alias	FOXO2; AF6q21; FKHL1; FOXO3A; FKHL1P2
Species	Human
Gene ID (Human)	2309
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 574-673 of human FOXO3A (NP_001446.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O43524
Protein Weight	71kDa
Shipping	Ice bag