

Anti-Foxo3A Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Wvh4

Item Number: CM0017291

Introduction

Rabbit polyclonal antibody against mouse FOXO3A, validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human and rat. Ideal for apoptosis, autophagy, and transcription research.

[Learn More](#)

Parameter	Value
Product Name	FOXO3A Rabbit pAb
Remarks/Alias	Fkhr2, FKHL1, Foxo3a, 1110048B16Rik, 2010203A17Rik, FOXO3A
Species	Mouse
GeneID (Human)	2309
GeneID	56484
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 558-635 of mouse FOXO3A (NP_062714.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9WVH4
Protein Weight	71kDa
Shipping	Ice bag