

Anti-Foxp1 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9H334

Item Number: CM0006847



Introduction

FOXP1 rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat FOXP1. Suitable for transcription regulation, B-cell development, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	FOXP1 Rabbit mAb
Remarks/Alias	MFH; QRF1; 12CC4; hFKH1B; HSPC215; FOXP1
Species	Human
GeneID	27086
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 550-677 of human FOXP1 (NP_116071.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H334
Protein Weight	75kDa
Shipping	Ice bag