

## Anti-Foxp2 Polyclonal Antibody For Wb, Elisa - 015409

Item Number: CM0030570



### Introduction

Rabbit polyclonal FOXP2 antibody validated for WB and ELISA in human, mouse, and rat samples. Targets the forkhead box protein P2 involved in speech and language development. Suitable for lung epithelium and neural tissue studies.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | FOXP2 Rabbit pAb                                                                                                   |
| Remarks/Alias    | SPCH1; CAGH44; TNRC10; FOXP2                                                                                       |
| Species          | Human                                                                                                              |
| GeneID (human)   | 93986                                                                                                              |
| GeneID           | 93986                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 590-715 of human FOXP2 (NP_055306.1) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | O15409                                                                                                             |
| Protein Weight   | 80kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |