

Anti-Foxp3 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Bzs1

Item Number: CM0004389



Introduction

High-quality FOXP3 Rabbit Monoclonal Antibody validated for Western blot, IHC-P, and ELISA applications. Reacts with human, mouse, and rat. Ideal for Treg cell research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	FOXP3 Rabbit mAb
Remarks/Alias	JM2; AIID; IPEX; PIDX; XPID; DIETER
Species	Human
GeneID (Human)	50943
GeneID	50943
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BZS1
Protein Weight	47kDa
Shipping	Ice bag