

Anti-Foxp3 Monoclonal Antibody For Wb, Ip, Elisa - Q9Bzs1

Item Number: CM0005107



Introduction

FOXP3 Rabbit mAb (Q9BZS1) validated for WB, IP, ELISA with human cross-reactivity; targets a forkhead transcription factor crucial for Treg development and immune homeostasis.

[Learn More](#)

Parameter	Value
Product Name	FOXP3 Rabbit mAb
Remarks/Alias	JM2; AIID; IPEX; PIDX; XPID; DIETER; FOXP3
Species	Human
GeneID (Human)	50943
GeneID	50943
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	Q9BZS1
Protein Weight	47kDa
Shipping	Ice bag