

## Anti-Foxp3 Polyclonal Antibody For Wb, Elisa - Q99Jb6

Item Number: CM0028464



### Introduction

Rabbit polyclonal antibody targeting mouse FOXP3, validated for Western blot and ELISA. FOXP3 is essential for regulatory T-cell development and immune homeostasis, making this antibody suitable for immunology and autoimmunity studies.

[Learn More](#)

| Parameter        | Value                                                                                               |
|------------------|-----------------------------------------------------------------------------------------------------|
| Product Name     | FOXP3 Rabbit pAb                                                                                    |
| Remarks/Alias    | sf; JM2; scurfin                                                                                    |
| Species          | Mouse                                                                                               |
| Gene ID (Human)  | 50943                                                                                               |
| Gene ID          | 20371                                                                                               |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse FOXP3 (NP_473380.1) |
| Source           | Rabbit                                                                                              |
| Category         | Polyclonal Antibodies                                                                               |
| Application      | WB, ELISA                                                                                           |
| Cross Reactivity | Mouse                                                                                               |
| SWISS            | Q99JB6                                                                                              |
| Protein Weight   | 47kDa                                                                                               |
| Shipping         | Ice bag                                                                                             |