

Anti-Foxp3 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Bzs1

Item Number: CM0018197



Introduction

Rabbit polyclonal antibody against human FOXP3, validated for Western blot and ELISA. Cross-reacts with rat. Ideal for regulatory T-cell biology, immune homeostasis, and IPEX disease research.

[Learn More](#)

Parameter	Value
Product Name	FOXP3 Rabbit pAb
Remarks/Alias	JM2; AIID; IPEX; PIDX; XPID; DIETER; FOXP3
Species	Human
GeneID (Human)	50943
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 70-180 of human FOXP3 (NP_054728.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q9BZS1
Protein Weight	47kDa
Shipping	Ice bag