

Anti-Foxq1 Rabbit Polyclonal Antibody For Ihc-P, Elisa - Q9C009

Item Number: CM0021224

Introduction

FOXQ1 Rabbit Polyclonal Antibody is designed for IHC-P and ELISA applications. It targets human FOXQ1, a transcription factor involved in hair follicle differentiation. Ideal for research in dermatology and developmental biology.

[Learn More](#)

Parameter	Value
Product Name	FOXQ1 Rabbit pAb
Remarks/Alias	FOXQ1; HFH1; forkhead box Q1
Species	Human
GeneID	94234
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FOXQ1 (NP_150285.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q9C009
Protein Weight	41kDa
Shipping	Ice bag