

Anti-Foxs1 Polyclonal Antibody For Wb, Elisa - O43638

Item Number: CM0012861

Introduction

Rabbit anti-FOXS1 polyclonal antibody for WB and ELISA. Validated in human and rat. FOXS1 is a transcriptional repressor involved in testicular vasculature organization. Ideal for research use.

[Learn More](#)

Parameter	Value
Product Name	FOXS1 Rabbit pAb
Remarks/Alias	FKHL18; FREAC10; FOXS1
Species	Human
GeneID (Human)	2307
GeneID	2307
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 121-330 of human FOXS1 (NP_004109.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	O43638
Protein Weight	35kDa
Shipping	Ice bag