


Anti-Fpr2 Rabbit Polyclonal Antibody For If/Icc And Elisa - P25090

Item Number: CM0022772



Introduction

Rabbit polyclonal antibody against human FPR2 (N-formyl peptide receptor 2) validated for IF/ICC and ELISA. Reacts with human, mouse, rat. Target involved in innate immunity and inflammation signaling.

[Learn More](#)

Parameter	Value
Product Name	FPR2 Rabbit pAb
Remarks/Alias	ALX; ALXR; HM63; FMLPX; FPR2A; FPRH1; FPRH2; FPRL1; LXA4R; FMLP-R-II; FPR2
Species	Human
Gene ID	2358
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human FPR2 (NP_001453.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25090
Protein Weight	39kDa
Shipping	Ice bag