

Anti-Fps/Fdps Monoclonal Antibody For Wb, Elisa - P14324

Item Number: CM0004634



Introduction

Rabbit monoclonal antibody targeting human FPS/FDPS, suitable for WB and ELISA. Reacts with human, mouse, rat. Farnesyl pyrophosphate synthase, key enzyme in isoprenoid biosynthesis, ~48kDa.

[Learn More](#)

Parameter	Value
Product Name	FPS/FDPS Rabbit mAb
Remarks/Alias	FPS; FPPS; POROK9; FPS/FDPS
Species	Human
GeneID (Human)	2224
GeneID	2224
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 67-419 of human FPS/FDPS. (NP_001995.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14324
Protein Weight	48kDa
Shipping	Ice bag