

Anti-Fra1 Monoclonal Antibody For Wb, Elisa - P15407

Item Number: CM0007722

Introduction

High-quality anti-FRA1 rabbit monoclonal antibody validated for Western blot (WB) and ELISA. Detects human Fos-related antigen 1 (FOSL1/FRA1), a ~29 kDa nuclear transcription factor. Suitable for cancer and cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	FRA1 Rabbit mAb
Remarks/Alias	FRA; FRA1; fra-1
Species	Human
GeneID (Human)	8061
GeneID	8061
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FRA1 (NP_005429.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P15407
Protein Weight	29 kDa
Shipping	Ice bag