


Anti-Frataxin (Fxn) Polyclonal Antibody For Wb - Q16595

Item Number: CM0018928



Introduction

Rabbit polyclonal antibody against FXN/Frataxin. Validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for studies on iron-sulfur cluster biogenesis and Friedreich ataxia.

[Learn More](#)

Parameter	Value
Product Name	FXN / Frataxin Rabbit pAb
Remarks/Alias	FA; X25; CyaY; FARR; FRDA; FXN / Frataxin
Species	Human
GeneID (Human)	2395
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 42-210 of human FXN / Frataxin (NP_000135.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16595
Protein Weight	23kDa
Shipping	Ice bag