


Anti-Frataxin/Fxn Rabbit Monoclonal Antibody [Kd Validated] For Wb, If/Icc, Ihc-P, Elisa - Q16595

Item Number: CM0005054



Introduction

Rabbit monoclonal antibody against human Frataxin (FXN), KD validated. Suitable for WB, IF/ICC, IHC-P, ELISA. Targets mitochondrial protein involved in iron-sulfur cluster assembly and iron homeostasis. For research use only.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] Frataxin/FXN Rabbit mAb
Remarks/Alias	FA; X25; CyaY; FARR; FRDA
Species	Human
Gene ID	2395 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 42-210 of human FXN/Frataxin (NP_000135.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q16595
Protein Weight	23kDa
Shipping	Ice bag