

Anti-Frem2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q5Szk8

Item Number: CM0019799

Introduction

Rabbit polyclonal antibody targeting human FREM2. Validated for WB, IF/ICC, ELISA. Cross-reacts with human and mouse. Ideal for studying epithelial integrity and cell adhesion. SWISS: Q5SZK8.

[Learn More](#)

Parameter	Value
Product Name	FREM2 Rabbit pAb
Remarks/Alias	CRYPTOP; FRASRS2; FREM2
Species	Human
Gene ID (Human)	341640
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 47-300 of human FREM2 (NP_997244.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q5SZK8
Protein Weight	351kDa
Shipping	Ice bag