

Anti-Frmd6 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q96Ne9

Item Number: CM0027846

Introduction

CamelBio offers FRMD6 rabbit polyclonal antibody (Q96NE9) for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat FRMD6, a FERM domain protein at membrane and cytoplasm.

[Learn More](#)

Parameter	Value
Product Name	FRMD6 Rabbit pAb
Remarks/Alias	EX1; Willin; C14orf31; c14_5320; FRMD6
Species	Human
Gene ID (Human)	122786
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 82-183 of human FRMD6 (NP_001253975.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96NE9
Protein Weight	72kDa
Shipping	Ice bag