

Anti-Frs2 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q8Wu20

Item Number: CM0005175



Introduction

FRS2 rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat FRS2 (SNT-1). Ideal for FGFR signaling pathway research and assay development.

[Learn More](#)

Parameter	Value
Product Name	FRS2 Rabbit mAb
Remarks/Alias	SNT; SNT1; FRS1A; FRS2A; SNT-1; FRS2alpha; FRS2
Species	Human
Gene ID	10818
Immunogen	Recombinant protein of human FRS2
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WU20
Protein Weight	57kDa
Shipping	Ice bag