

Anti-Frs2 Rabbit Polymab® For Wb, Ihc-P, If/Icc, Elisa - Q8Wu20

Item Number: CM0005505



Introduction

FRS2 Rabbit PolymAb® is a rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, ELISA with reactivity to human, mouse, rat. FRS2 is an adaptor protein linking FGFR and NGF receptors to MAP kinase and PI3K signaling pathways. Ideal for signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	FRS2 Rabbit PolymAb®
Remarks/Alias	SNT; SNT1; FRS1A; FRS2A; SNT-1; FRS2alpha
Species	Human
Gene ID	10818
Immunogen	Recombinant protein (recombinant protein of human FRS2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WU20
Protein Weight	57 kDa
Shipping	Ice bag