

Anti-Frzb Rabbit Polyclonal Antibody For Wb, Elisa - Q92765

Item Number: CM0030435

Introduction

Rabbit polyclonal antibody against human FRZB (SFRP3, FRZB-1). Validated for WB and ELISA, cross-reacts with rat and mouse. Targets a secreted Wnt signaling modulator involved in skeletogenesis.

[Learn More](#)

Parameter	Value
Product Name	FRZB Rabbit pAb
Remarks/Alias	FRE, OS1, FZRB, hFIZ, FRITZ, FRP-3, FRZB1, SFRP3, SFRP3, FRZB-1, FRZB-PEN, FRZB
Species	Human
GeneID (human)	2487
GeneID	2487
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 166-325 of human FRZB (NP_001454.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat, Mouse
SWISS	Q92765
Protein Weight	36kDa
Shipping	Ice bag