

Anti-Fsh Alpha/Beta Rabbit Monoclonal Antibody For If-P, Ihc-P, Elisa - P01215 / P01225

Item Number: CM0004357



Introduction

Rabbit monoclonal antibody targeting FSH alpha/beta chains, validated for IF-P, IHC-P, and ELISA applications. Specific for human samples. Ideal for reproductive hormone research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	FSH alpha/beta Rabbit mAb
Remarks/Alias	HCG; LHA; FSHA; GPA1; GPHa; TSHA; GPHA1; CG-ALPHA; FSHB
Species	Human
GeneID (Human)	1081/2488
GeneID	1081/2488
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P01215/P01225
Protein Weight	13kDa/15kDa
Shipping	Ice bag