

Anti-Fshb Polyclonal Antibody For Wb, Elisa - P01225

Item Number: CM0016695



Introduction

Rabbit polyclonal antibody against human FSHB for WB and ELISA applications. Cross-reacts with mouse and rat. Suitable for reproductive biology research and diagnostic assay development.

[Learn More](#)

Parameter	Specification
Product Name	FSHB Rabbit pAb
Remarks/Alias	HH24; FSHB
Species	Human
GeneID (Human)	2488
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-129 of human FSHB (NP_000501.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P01225
Protein Weight	15kDa
Shipping	Ice bag