

Anti-Fsp1/S100A4 Rabbit Polyclonal Antibody For Wb - P26447

Item Number: CM0020399



Introduction

Anti-FSP1/S100A4 Rabbit pAb validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for studying S100A4 in cell motility, metastasis, and apoptosis. High-quality IVD raw material.

[Learn More](#)

Parameter	Value
Product Name	FSP1/S100A4 Rabbit pAb
Remarks/Alias	42A; 18A2; CAPL; FSP1; MTS1; P9KA; PEL98; FSP1/S100A4
Species	Human
Gene ID	6275
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-101 of human FSP1/S100A4 (NP_002952.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26447
Protein Weight	12kDa
Shipping	Ice bag