

Anti-Fstl4 Polyclonal Antibody For Wb, If/Icc, Elisa - Q6Mzw2

Item Number: CM0020537

Introduction

FSTL4 Rabbit pAb is a polyclonal antibody validated for Western blot, IF/ICC, and ELISA. It detects human, mouse, and rat FSTL4, a secreted follistatin family protein. Ideal for extracellular signaling and developmental biology studies.

[Learn More](#)

Parameter	Value
Product Name	FSTL4 Rabbit pAb
Remarks/Alias	FSTL4
Species	Human
GeneID (human)	23105
GeneID	23105
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-250 of human FSTL4 (NP_055897.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6MZW2
Protein Weight	93kDa
Shipping	Ice bag