

Anti-Fstl5 Polyclonal Antibody For Wb, If-P, Elisa - Q8N475

Item Number: CM0020649

Introduction

High-quality rabbit polyclonal antibody against human FSTL5, validated for WB, IF-P, and ELISA. Cross-reacts with mouse and rat. Ideal for studying follistatin-related protein 5 in various samples.

[Learn More](#)

Parameter	Value
Product Name	FSTL5 Rabbit pAb
Remarks/Alias	FSTL5
Species	Human
Gene ID (Human)	56884
Gene ID	56884
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-250 of human FSTL5 (NP_064501.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N475
Protein Weight	96kDa
Shipping	Ice bag